Penumbra Embolization Coils MRI Safety Information: MR Conditional

Non-clinical testing has demonstrated the coil is MR Conditional. Coil arrays consisting of one or more coil implants can be scanned safely immediately after implantation under the following conditions.

Parameter	Condition
Static Magnetic Field Strength	1.5T or 3T
Maximum Spatial Field Gradient	3000 Gauss/cm (30 T/m) or less
RF Coil Type	Any
Maximum Whole Body SAR	2.0 W/kg
Maximum Head SAR	3.2 W/kg

RF Heating Information

Under the scanning conditions defined above, the coils are expected to produce an estimated in-vivo temperature rise < 3°C during 15 minutes of continuous MR scanning. This temperature rise is conservative, as it does not include the cooling effects of blood flow.

Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the Penumbra coil arrays. Therefore, it may be necessary to optimize parameters for the presence of this implant. In non-clinical testing:

Penumbra Coil 400™, POD400™, Ruby™ Coil, POD™ and LP Coil:

- Artifact associated with a Penumbra coil array was assessed using spin echo and gradient echo sequences. The maximum artifact distance beyond the implant was 3 mm for the spin echo sequence and 8 mm for the gradient echo sequence.
- Artifact associated with a Penumbra coil array in the form of an 8 mm diameter sphere was assessed using a clinical MRA sequence. The maximum artifact distance beyond the implant was 2 mm using this sequence.

SMART COIL™:

- The distance of maximum extent of the artifact beyond the metal of the array was 3 mm for the spin echo sequence and 18 mm for the gradient echo sequence.

Patient's Name

Treating Physician

Phone Number

Hospital information

Date of treatment or other information

If you have any questions about MRI scans or Penumbra Embolization Coils, contact your physician.

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