









Indigo™ System

Coronary Mechanical Thrombectomy

CAT™ RX

Coronary Data Overview

	Initial Experience	Stony Brook Experience	UH Cleveland Experience ^e	CHEETAH	
	72	34	123	400	Patients
	Retrospective Multicenter Investigator Initiated	Retrospective Single-center Investigator Initiated Stony Brook University, NY	Retrospective Single-center Investigator Initiated University Hospitals Cleveland, OH	Prospective Multicenter Frontline use of CAT RX Penumbra Sponsored	Study Design
	RCA, LAD, LCx	RCA, LAD, LCx, SVG, Ramus	RCA, LAD, LCx, SVG, LM, OM	RCA, LAD, LCx	Primary Vessels Treated
	86.1%	88%	—	85% (TIMI 2–3 flow)	Post-CAT RX TIMI 3 Flow
	97.2%	—	90.2% 70% of patients who did not achieve TIMI 3 flow received balloon dilatation prior to CAT RX	97.5%	Final TIMI 3 Flow
	0^a	0^b	0^d	0^e device-related	Ischemic Stroke
	<p>In our initial experience, aspirating thrombus from ACS patients using the Indigo CAT RX Aspiration System before PCI was safe and effective for reducing thrombus burden and restoring flow.</p>	<p>Penumbra CAT RX system appears to be safe and highly effective at thrombus removal in the acute setting without increased stroke risk seen with manual aspiration.</p>	<p>We strongly believe that the Penumbra CAT RX device should be used before any balloon angioplasty in patients with high thrombus burden...</p>	<p>Sustained power aspiration with CAT RX is safe for the removal of high thrombus burden with low rates of distal embolization and improved myocardial perfusion.</p>	Key Takeaways
	<p>a. There was 1 (1.4%) incident of hemorrhagic stroke reported. Timeframe: 30 days post-procedure. Mathews SJ, et al. Initial experience with a continuous mechanical aspiration system for thrombus removal before percutaneous coronary intervention. <i>Catheter Cardiovasc Interv.</i> 2022 Nov;100(6):950–954. doi:10.1002/ccd.30389. Epub 2022 Oct 2.</p>	<p>b. Timeframe: post-procedure. Gilchrist Jr IC, Fordham MJ, Pyo R, et al. Mechanical aspiration thrombectomy using the penumbra CAT RX system for patients presenting with acute coronary syndrome. <i>Cardiovasc Revasc Med.</i> 2021 Jul 2;S1553–8389(21)00501–7. doi:10.1016/j.carrev.2021.06.130.</p>	<p>c. 11 patients (8.9%) experienced an adverse cardiovascular event (3 MI, 8 heart failure) and there were 14 patient (11.3%) deaths within 6 months. d. Timeframe: not specified. Tashish N, Chami T, Dong T, et al. Routine use of the “Penumbra” thrombectomy device in myocardial infarction: a real-world experience—ROPUST study. <i>J Interv Cardiol.</i> 2022 Mar 26;2022:5692964. doi:10.1155/2022/5692964.</p>	<p>e. There were 3 (7.7%) incidence of non-device related strokes as adjudicated by Independent Medical Reviewer. Timeframe: 30 days post-procedure. Mathews SJ, Parkh SA, Wu W, et al. Sustained mechanical aspiration thrombectomy for high thrombus burden coronary vessel occlusion: the multicenter CHEETAH study. <i>Circ Cardiovasc Interv.</i> 2023 Feb;16(2):e012433. doi:10.1161/CIRCINTERVENTIONS.122.012433. Epub 2023 Feb 21. The safety and effectiveness of this device for use in the treatment of ST-Elevation Myocardial Infarction (STEMI) has not been established. Complications from the use of this device in this manner could lead to death, permanent impairment, and/or the need for emergency medical intervention. Penumbra ENGINE™ and Penumbra Pump MAX™ were both used during CHEETAH Study.</p>	



Maximized Clot Removal

Large .044" (1.12 mm) Lumen




Penumbra ENGINE™

Rapid Exchange

.014" (.356 mm) Wire Compatibility


Advanced Trackability

with Neuro-Tracking Technology




Enhanced Deliverability

Proximal Laser-Cut Hypotube



Sustained Power Aspiration vs. Manual Syringe-Based Aspiration



112%

more aspiration efficiency compared with a syringe^h



Watch the Video Now to see Sustained Power Aspiration with Penumbra ENGINE

f. Tests performed and data on file at Penumbra, Inc. Bench test results may not be indicative of clinical performance.
 g. Tests performed and data on file at Penumbra, Inc. Bench test results may not be indicative of clinical performance. Data obtained using a 30 cc syringe connected to manual aspiration catheter with vacuum gauge attached to proximal barrel of syringe and continuous filling of syringe.
 h. Aspiration efficiency defined as volume under graph.



Penumbra



penumbrainc.com

Penumbra, Inc. USA

One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

Penumbra Europe GmbH

Am Borsigturm 44
13507 Berlin
Germany
T +49 30 2005 676-0
F +49 30 2005 676-10
de-order@penumbrainc.com
de-info@penumbrainc.com

Penumbra Neuro Australia Pty Ltd

55 Kirby Street
Rydalmere NSW 2116
Australia
T +61-1300 817 025
F +61-1300 817 026
order.anz@penumbrainc.com

**Penumbra Latin America
Distribuidora de Equipamentos
e Produtos Médicos Ltda**

Av. Brigadeiro Faria Lima 1336
Cj 82, CEP 01451-001
São Paulo/SP, Brazil
T +55 11 2883-5825
order.la@penumbrainc.com



For the complete
Penumbra IFU Summary
Statements, please
scan QR code or visit:
peninc.info/risk

Background image used under license from shutterstock.com. The clinical results presented herein are for informational purposes only, and may not be predictive for all patients. Individual results may vary depending on patient-specific attributes and other factors. Tests performed and data on file at Penumbra, Inc. Bench test results may not be indicative of clinical performance.

Product availability varies by country. Caution: Federal (USA) law restricts these devices to sale by or on the order of a physician. Prior to use, please refer to the Instructions for Use for complete product indications, contraindications, warnings, precautions, potential adverse events, and detailed instructions for use. Please contact your local Penumbra representative for more information.

Copyright ©2023 Penumbra, Inc. All rights reserved. The Penumbra P logos, CAT, Indigo, Penumbra ENGINE, and MAX are registered trademarks or trademarks of Penumbra, Inc. in the USA and other countries. 25562, Rev. C 12/23 UNI