

Results of MISTIE III

Presented at International Stroke Conference
February 7, 2019; Honolulu, HI, USA

**Evaluating Image-Guided, Minimally Invasive
Surgery for ICH: MISTIE III Results**

Presented by Daniel F. Hanley, MD

**MISTIE III Surgical Results: Efficiency of
Hemorrhage Removal Determines mRS**

Presented by Issam A. Awad, MD

Publication of Study Results

**Efficacy and safety of minimally invasive surgery
with thrombolysis in intracerebral haemorrhage
evacuation (MISTIE III): a randomised, controlled,
open-label, blinded endpoint phase 3 trial.**

Hanley DF, Thompson RE, Rosenblum M, Yenokyan G, Lane K, McBee N, Mayo SW, Bistran-Hall AJ, Gandhi D, Mould WA, Ullman N, Ali H, Carhuapoma JR, Kase CS, Lees KR, Dawson J, Wilson A, Betz JF, Sugar EA, Hao Y, Avadhani R, Caron JL, Harrigan MR, Carlson AP, Bulters D, LeDoux D, Huang J, Cobb C, Gupta G, Kitagawa R, Chicoine MR, Patel H, Dodd R, Camarata PJ, Wolfe S, Stadnik A, Money PL, Mitchell P, Sarabia R, Harnof S, Barzo P, Unterberg A, Teitelbaum JS, Wang W, Anderson CS, Mendelow AD, Gregson B, Janis S, Vespa P, Ziai W, Zuccarello M, Awad IA; MISTIE III Investigators. *Lancet*. 2019 Feb 6. pii: S0140-6736(19)30195-3. doi: 10.1016/S0140-6736(19)30195-3.

Removing ICH with MIS is Safe

Safety Outcomes ¹			
	MISTIE (n=255)	Standard medical care (n=251)	p value
Died within 0-7 days	2 (1%)	10 (4%)	.018
Died within 0-30 days	24 (9%)	37 (15%)	.066
Died within 0-180 days	39 (15%)	57 (23%)	.033

69%

Average clot reduction


59%

Patients with ≤ 15 mL EOT

145/246 patients randomized to MISTIE procedure

Minimally Invasive Surgical Reduction of Clot Volume to ≤ 15 mL Produces Better Functional Outcomes

Dichotomized Outcome with ≤ 15 mL EOT (Residual) Volume at Day 365 (as Treated)



- Cases with ≤ 15 mL EOT (residual) volume had lower mortality and better functional outcome
- Reduction beyond 15 mL significantly increased the chance of a good functional outcome, by 10% for each additional mL hematoma removed (p = .002)
- Further volume reduction beyond 70% removal carries a significant benefit, with 6% improvement in chance of achieving mRS 0-3 per additional mL removed (p = .002)

1. Frequency of survival is modestly improved without a "price" in surgical risk or vegetative state.

MISTIE III Trial Design

Study Protocol

Inclusion Criteria

- ICH \geq 30 mL
- ICH/IVH/IVH catheter tract/BP stability
- Randomize 12 to 72 hours post onset
- Age \geq 18 years
- Historical modified Rankin Scale score \leq 1

Exclusion Criteria

- Vascular defect R/O by CTA
- Infratentorial hemorrhage; evidence of brain stem involvement; large IVH
- Anticoagulation required; irreversible platelet count $<$ 100,000 or INR $>$ 1.4
- Uncontrollable systemic bleeding
- Other comorbidity preventing use of thrombolytic therapy or follow-up

Baseline Characteristics (mITT Group)

	MISTIE (n=250)	Standard medical care (n=249)
Sex		
Men	159 (64%)	146 (59%)
Women	91 (36%)	103 (41%)
GCS score at randomization*		
3–8	64 (26%)	63 (25%)
9–12	111 (44%)	108 (43%)
13–15	75 (30%)	78 (31%)
Diagnostic CT at presentation		
ICH volume (mL)	42.7 (30.4–54.5)	41.5 (30.9–55.3)
Stability CT (last CT before randomization)		
ICH volume (mL)	45.8 (35.4–59.6)	45.3 (35.4–57.2)
Clot location		
Deep	163 (65%)	144 (58%)
Lobar	87 (35%)	105 (42%)
mRS score before stroke		
0	230 (92%)	233 (94%)
1	20 (8%)	16 (6%)

Penumbra, Inc. USA
One Penumbra Place
Alameda, CA 94502
USA
1.888.272.4606
T 1.510.748.3200
F 1.510.748.3232
order@penumbrainc.com
info@penumbrainc.com

Penumbra Europe GmbH
Am Borsigturm 44
13507 Berlin
Germany
T +49 30 2005 676-0
T +49 30 2005 676-10
de-order@penumbrainc.com
info@penumbrainc.de

Penumbra Neuro Australia Pty Ltd
Suite 3, Level 5, 1 Oxford Street
Darlinghurst NSW 2010
Australia
T +61-1300 817 025
F +61-1300 817 026
order.anz@penumbrainc.com

The clinical results presented herein are for informational purposes only, and may not be predictive for all patients. Individual results may vary depending on a variety of patient-specific attributes. Please contact your local Penumbra representative for more information.

Copyright ©2019 Penumbra, Inc. All rights reserved. The Penumbra P logo is a registered trademark or trademark of Penumbra, Inc. in the USA and other countries. 15639, Rev. A 03/19 UNI